SNAP* Cortisol Test Kit

For use only with the IDEXX SNAP* Reader Analyzer and the IDEXX SNAPshot Dx* Analyzer.

For veterinary use only.

English Version

Overview

The SNAP Cortisol Test Kit is an enzyme-linked immunosorbent assay (ELISA) for the quantitative measurement of cortisol in dog serum. For in vitro use only.

SNAP Cortisol Test Ranges

Dynamic ranges:

- 2.5–30 µg/dL (69–828 nmol/L) for post-ACTH stimulation, Cushing’s-suspected dogs and therapeutic monitoring > 10 µg/dL (>276 nmol/L)
- 0.5–10 µg/dL (14–276 nmol/L) for all other tests

Precautions and Warnings


- Run the SNAP calibration device at the beginning of each day (SNAP Reader analyzer only).
- This test cannot be visually interpreted. The IDEXX SNAP Reader analyzer or SNAPshot Dx analyzer must read all tests.
- Use a separate pipette tip and sample tube for each test.
- Do not use components after their expiration dates.
- Do not mix components from different lots of test kits.
- Do not use a device that does not have two blue stripes visible before adding the sample.
- Keep the device in its sealed foil pack at room temperature (18–27°C/64–80°F) for at least 10 minutes before use.
- Prewarming of the conjugate solution is not required.
- Keep the SNAP device horizontal during the entire test procedure.**
- IMPORTANT: Do not depress the activator until instructed.**
- Insert the device into the analyzer immediately after activation.
- WARNING:** Conjugate – H316/P332+P313/EUH208. Causes mild skin irritation. If skin irritation occurs. Get medical advice/attention. Contains Kathon. May produce an allergic reaction.

Storage

- SNAP devices and test reagents must be stored at 2–7°C (36–45°F).
- Do not heat or microwave the kit components.
- Take out only the kit components you will need.
- Do not leave kit components at room temperature for more than 8 consecutive hours.

Kit Components

- Conjugate solution
- SNAP test devices
- Disposable pipette tips for 300-µL manual pipettor
- Disposable pipette tips for the 25-µL and 100-µL manual pipettors (for use with the SNAPshot Dx analyzer only)
- Disposable sample tubes
- Workstation

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Componentes Required (Not Provided in Kit)

SNAP Reader analyzer users:

- VetTest* pipette tips
- 300-µL manual pipettor (provided with the SNAP Reader analyzer)

SNAPshot Dx analyzer users:

The following required components are provided with the SNAPshot Dx analyzer:

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Sample Information

- Do not use plasma or whole blood. Use serum, either fresh, previously frozen or stored at 2–7°C (36–45°F) for up to 7 days. To store longer than 7 days, freeze the sample (-20°C/4°F or below). Remove serum from the red blood cells before freezing. Bring previously refrigerated or frozen samples to room temperature 18–27°C (64–80°F) before running assays.
- SNAP cortisol sample size requirements:
 - 25 µL for post-ACTH stimulation, Cushing’s-suspected dogs and therapeutic monitoring > 10 µg/dL (>276 nmol/L)
 - 100 µL for all other tests
- Centrifuge previously frozen or older samples before use.
- Results of this test are not affected by hemolysis or lipemia.

Test Procedure

SNAP Reader analyzer users:

The VetTest analyzer provides on-screen instructions throughout the analysis. For complete instructions, refer to the Running Quantitative SNAP Tests chapter in the SNAP Reader Operator’s Manual.

SNAPshot Dx analyzer users:

The SNAPshot Dx analyzer provides on-screen instructions throughout the analysis. For complete instructions, refer to the Using the SNAPshot Dx Analyzer chapter in the SNAPshot Dx Operator’s Manual.

Test	Sample Volume	Conjugate Volume
Baseline Cortisol	100 µL	300 µL
Post-Dexamethasone Cortisol (post low- or high-dose dexamethasone)	100 µL	300 µL
Post-ACTH Stim Cortisol: Cushing’s Suspected	25 µL	300 µL
Post-ACTH Stim Cortisol: Addison’s Suspected	100 µL	300 µL
Post-ACTH Stim Cortisol: Therapeutic Monitoring	100 µL	300 µL
Post-ACTH Stim Cortisol: Therapeutic Monitoring >10 µg/dL	25 µL	300 µL

Interpreting Test Results

The recommendations below are intended to provide general guidance only. As with any diagnosis, you should use clinical discretion with each patient based on a complete evaluation of the patient, including history, physical examination and complete laboratory profile.

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A. ACTH Stimulation Test

Pre-ACTH	Interpretation	Post-ACTH	Interpretation
<2 µg/dL	Must evaluate in conjunction with post-ACTH result. If both results are <2 µg/dL, results are consistent with hypoadrenocorticism. Begin treatment with mineralocorticoid and/or glucocorticoid as appropriate.	<2 µg/dL	Ideally, should be evaluated in conjunction with pre-ACTH result. If both results are <2 µg/dL, results are consistent with hypoadrenocorticism. Begin treatment with mineralocorticoid and/or glucocorticoid as appropriate.
2–6 µg/dL	Normal	2–6 µg/dL	Inconclusive
>22 µg/dL	Consistent with Cushing’s syndrome; perform high-dose dexamethasone suppression.	6–18 µg/dL	Normal
		18–22 µg/dL	Equivocal; Cushing’s syndrome possible
		>22 µg/dL	Consistent with Cushing’s syndrome; perform high-dose dexamethasone suppression to discriminate between PDH and ATH, ACTH level and/or abdominal ultrasound.

B. Low-Dose Dexamethasone Suppression Test

4-hour cortisol level	8-hour cortisol level	Interpretation
Not necessary	<1 µg/dL	Normal
1–1.5 µg/dL	1–1.5 µg/dL	Inconclusive, consider repeating in 8–12 weeks
>1.5 µg/dL and >50% of baseline	>1.5 µg/dL and >50% of baseline	Consistent with Cushing’s syndrome; perform high-dose dexamethasone suppression, endogenous ACTH concentration and/or abdominal ultrasound to discriminate between PDH and ATH.
<1.5 µg/dL or <50% of baseline	>1.5 µg/dL and >50% of baseline	Consistent with PDH
>1.5 µg/dL or >50% of baseline	>1.5 µg/dL and <50% of baseline	Consistent with PDH
<1.5 µg/dL or <50% of baseline	>1.5 µg/dL and <50% of baseline	Consistent with PDH

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


C. High-Dose Dexamethasone Suppression Test

4-hour cortisol level	8-hour cortisol level	Interpretation
<1.5 µg/dL or <50% of baseline	>1.5 µg/dL and >50% of baseline	Consistent with PDH
>1.5 µg/dL and >50% of baseline	<1.5 µg/dL or <50% of baseline	Consistent with PDH
<1.5 µg/dL or <50% of baseline	<1.5 µg/dL or <50% of baseline	Consistent with PDH
>1.5 µg/dL and >50% of baseline	>1.5 µg/dL and >50% of baseline	Additional testing required to differentiate PDH from ATH; consider measuring plasma ACTH levels and/or performing an abdominal ultrasound.

Symbol Descriptions		
 Use by date		Authorized Representative in the European Community
 Batch Code (Lot)		Consult instructions for use
 Catalog Number		In vitro diagnostic
 Temperature limitation		
 Manufacturer		

IDEXX Customer Support

US/Canada/Latin America: 1-800-248-2483

Europe: 00800 1234 3399

New Zealand: 0800-102-084

Australia: 1300 44 33 99

Asia: 0800-291-018

If possible, please have your kit lot number available.

*Valid only in countries where mitotane is approved for use; trilostane may be used as an alternative where approved for use, in accordance with applicable protocols.

*SNAP, SNAPshot Dx and VetTest are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries. All other products and company names are trademarks of their respective holders. Patent information: idexx.com/patents.

© 2015 IDEXX Laboratories, Inc. All rights reserved • 06-03312-10

	IDEXX Laboratories, Inc. One IDEXX Drive, Westbrook, Maine 04092 USA	
	IDEXX Europe B.V. PO. Box 1334, NL-2130 EK Hoofddorp	

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© 2015 IDEXX Laboratories, Inc. Tutti i diritti sono riservati.