## Canine Leishmania Antibody Test Kit

For veterinary use only.

### SNAP® Leishmania

The Canine Leishmania Antibody Test Kit is an Enzyme-linked Immunosorbent assay for *in vitro* detection of *Leishmania donovani* or *Leishmania infantum* (*L. donovani* or *L. infantum*) antibodies in canine whole blood, serum, or plasma. Leishmaniasis, caused by parasites of the genus Leishmania, has wide distribution throughout the Mediterranean Basin. The parasites are transmitted by sand flies to many mammals including humans and canines. Dogs are likely an important reservoir for both the human and canine diseases. Leishmania primarily infect white blood cells leading to impaired immune function, blood disorders, and various clinical signs and skin lesions. A viscerale leishmaniasis include weight loss, fever, and various clinical signs and skin lesions. In dogs, diagnosis of visceral leishmaniasis can be made following observation of typical clinical signs and the measurement of a significant antibody titer to Leishmania. The Canine Leishmania Antibody Test employs purified Leishmania antigen to detect Leishmania-specific antibody. The SNAP test was optimized with the widely used immunofluorescent assay (IFA) method as the reference standard. A positive result on the SNAP test correlates to a significant IFA titer.<sup>1</sup> The SNAP test is an enzyme immunoassay technology developed by IDEXX Laboratories, Inc. In operation, the conjugate and test sample are mixed and added to the SNAP device. The device is then activated, releasing reagents stored within the device. Color development is proportional to the concentration of Leishmania-specific (*L. donovani* or *L. infantum*) antibodies in canine whole blood, serum or plasma. A Positive Control Spot develops color indicating assay reagents are active. Color development in the Negative Control Spot indicates reactivity unrelated to the disease or an improperly performed test.

### Precautions and Warnings

- Use a separate sample tube and transfer pipet for each test.
- During manufacturing the bioactive spots are dyed for quality control purposes. This does not interfere with the test results or interpretation.
- All waste should be properly decontaminated prior to disposal.
- Do not use components past expiration date.
- Do not mix components from kits with different lot numbers.
- The SNAP device must be in a horizontal position on a flat surface while performing the test.
- Do not use a SNAP device that has been activated prior to the addition of sample.

### Storage

- SNAP devices and test reagents must be stored at 2°C–8°C. All components must be at room temperature (15°C–25°C) before running the test—DO NOT HEAT.

### Kit Components

Item	Reagent	Amount	Amount
1.	<i>Leishmania</i> Antigen:HRPO Conjugate	1 x 6.0 mL	2 x 6.0 mL
2.	SNAP Device	10	30
Reagents contained in each device:			
Wash Solution	0.4 mL	0.4 mL	
Substrate Solution	0.6 mL	0.6 mL	
Other Components	transfer pipets, sample tubes, reagent rack		

### Sample Information

- Samples must be at room temperature (15°C–25°C) before beginning test procedure.
- Serum or plasma, either fresh, recently frozen or stored at 2°C–8°C, may be used in this test.
- Serum or plasma may be stored for up to 7 days at 2°C–8°C. For longer storage, sample should be frozen (-20°C or colder).
- Previously frozen or older samples must be centrifuged before use.
- Hemolyzed or lipemic samples will not affect results.
- EDTA or heparin in plasma will not affect results.
- Whole blood may be used. Whole blood must be anticoagulated (eg. EDTA, heparin) and may be used either fresh or after refrigeration (2°C–8°C) for up to one week.

### Test Procedure

- Using the pipet that is provided, transfer 2 drops of sample (whole blood, serum, or plasma) into sample tube (blue cap).
  - Holding bottle vertical, add 6 drops of conjugate to the sample tube (blue cap).
  - Cap the sample tube and mix thoroughly by inverting tube 3–5 times.
  - Place the device on a flat surface. Add contents of sample tube to Well, being careful not to splash contents outside of Sample Well. Sample will flow across Result Window, reaching Activate Circle in approximately 30–60 seconds. Some sample may remain in Sample Well.
- Watch carefully for sample or blue color in the Activate Circle.**
- When color **FIRST** appears in Activate Circle push Activator firmly until it is flush with the device body.
- 

**NOTE:** Some samples may not flow to the activate circle within 60 seconds, and, therefore, the circle may not turn color. In this case, press the Activator after sample has flowed across Result Window.

Keep the device horizontal to ensure accurate results.

Read test result at 6 minutes.

**NOTE:** positive control may develop sooner, but results are not complete until 6 minutes.

### Interpreting Test Result

To determine test results, read the reaction spots in the Result Window. Color development in Leishmania sample spot is proportional to the Leishmania-specific antibody concentration in the sample. If no color develops in the positive control spot, repeat the test.

### Negative Result

Only positive control spot develops color.

### Positive Result

Positive control spot and sample spot develop color.

### Reaction with Negative Control

The negative control spot serves as a safeguard against false positives and helps indicate that the assay has been run properly.

### Positive Result

If color in the sample spot is darker than negative control spot, result is positive for that sample spot.

### Invalid Results

1. **Background**  
If the sample is allowed to flow past the Activate Circle, background color may result. Some background color is normal. However, if colored background obscures test result, repeat the test.

2. **No Color Development**  
If positive control does not develop color, repeat the test.

### IDEXX Customer Support

USA/Canada 1-800-248-2483 • Europe 00-800-1234-3399 • Australia 1800 655 978

[idexx.com](http://idexx.com)

1. Georgi, JR, Georgi, ME. *Canine Clinical Parasitology*. Philadelphia, Lea and Febiger, 1992.

2. The SNAP test gave the same result as the IFA in 98% of the samples tested. IFA developed by Dr. Tesouro (Univ of Leon), study performed at Univ of Madrid, 108 IFA Positive Samples, 124 IFA Negative Samples.

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© 2009 IDEXX Laboratories, Inc. Alle Rechte vorbehalten.